MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007



TO10 & TO12

Threshold Contrast Phantoms

TO₁₀

These test objects are designed for quick quantitative assessments of image quality. The results can be plotted on a Threshold Detection Index Curve.

TO 10 is intended for use with Fluoroscopy systems.

108 details (12 sizes x 9 contrasts)
Size range 11mm to 0.25mm
Contrast range 0.012 to 0.930 @ 70kV, 1.0mm Cu filtration



TO12

These test objects are designed for quick quantitative assessments of image quality. The results can be plotted on a Threshold Detection Index Curve.

TO 12 is intended for use with Digital Spot Imaging Systems.

108 details (12 sizes x 9 contrasts)
Size range 11mm to 0.25mm
Contrast range 0.0043 to 0.954 @ 70kV, 1.0mm Cu filtration



MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007

TO10 & TO12

Threshold Contrast Phantoms

Product Radiograph




fig. 1 TO10 X-ray image

www.leedstestobjects.com

MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007

TO10 & TO12

Threshold Contrast Phantoms

Product Radiograph




fig. 2 TO12 X-ray image

www.leedstestobjects.com