

MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007



CR/DDR Set

Computed and Digital Radiography

A set of test objects designed to be used quickly and easily on an annual / commissioning basis to provide an ongoing check of imaging performance in computed radiography (CR) and digital radiography (DR/DDR), particularly those aspects which are most liable to deterioration.

The kit consists of:

- Resolution test pattern (type 18-010, up to 5LP/mm)
- TO20 threshold contrast test object (Threshold-contrast detail-detectability: 144 details (12 sizes x 12 contrasts), size range 11mm to 0.25mm, contrast range 0.0014 to 0.924 @ 75kV, 1.5mm Cu filtration)
- **TO M1** Geometry test object (Horizontal lines at 10mm pitch, increments notated along central axis)
- **TO MS4** Mesh Test Object for uniformity, blurring and stitching artefacts checks (nominal aperture size 0.50mm, nominal wire diameter 0.31mm)
- Small Lead block for image retention and erasure cycle efficiency checks
- Steel ruler
- Tape measure
- Adhesive tape
- 1.5mm Copper filter



www.leedstestobjects.com

tel: +44(0)1423 320 007

8.

CR/DDR Set

X-Ray Images



fig. 1 TO20 X-ray Image




fig. 2 TO M1 X-ray Image

www.leedstestobjects.com